



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
6-Methoxy-2-Naphthylacetic acid <i>aka Nabumetone metabolite</i>	Human Plasma	0.5-100 µg/mL	HPLC/FL	0.05 mL	EDTA	-20	Richmond	Inactive	LC184
3-Hydroxymandelic Acid	Human Urine	100-5000 ng/mL	LC/MS/MS	0.05 mL	None	-70	Madison	Active	P957.00
Acetaminophen  Acetaminophen (Acetaminophen sulfate) (p-Acetamidophenyl-β-D-glucuronide)	Human Plasma	0.025-10 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Active	LC12.1
	Human Urine	0.5-250 µg/mL	HPLC/UV	0.2 mL	None	-20	Richmond	Inactive	LC12.3
	Human Plasma	0.2-50 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Active	P089.03
	Human Plasma	0.2-50 µg/mL	HPLC/UV	0.1 mL	EDTA	-20	Madison	Active	P829.00
	Dog Plasma	0.2-100 µg/mL	HPLC/UV	0.3 mL	Sodium Heparin	-20	Madison	Active	P363.00
	Human Urine	0.2-200 µg/mL 0.5-500 µg/mL 1-1,000 µg/mL	LC/MS/MS	0.025 mL	None	-20	Madison	Active	P851.00
	Human Plasma	0.1-50 µg/mL 0.1-50 µg/mL 0.5-100 µg/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P852.00
Acyclovir	Human Plasma	5-1,000 ng/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC189
Adiponectin	Human Serum	185-23,695 ng/mL	ELISA	0.01 mL	None	-70	Richmond	Active	ICD287
Albendazole (Sulfoxide & Sulfone Metabolites)	Bovine Plasma	0.05-5 µg/mL 0.05-5 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Active	P562.02

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Albuterol	Human Plasma	0.20-20 ng/mL	LC/MS/MS	1.0 mL	Sodium Heparin	-20	Madison	Active	P478.01
	Human Plasma	2-4,000 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P951.00
	Human Urine	20-2,000 ng/mL	LC/MS/MS	0.1 mL	None	-20	Madison	Active	P465.01
	Dog Plasma	0.2-20 ng/mL	LC/MS/MS	1.0 mL	Sodium Heparin	-20	Madison	Active	P464.01
	Rat Plasma	0.4-400 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Madison	Inactive	P463.00
Alprazolam	Human Plasma	0.25-40 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Richmond	Active	LCMS225
Alprazolam ( $\alpha$ -Hydroxyalprazolam) (4-Hydroxyalprazolam)	Human Plasma	0.2-40 ng/mL 0.2-40 ng/mL 0.2-40 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Madison	Active	P335.06
Ambrisentan (4-Hydroxymethyl Ambrisentan)	Human Plasma	1-400 ng/mL 0.25 -100 ng/mL	LC/MS/MS	0.1 mL	Sodium Citrate	-20	Richmond	Active	LCMSC498
Amitriptyline (Nortriptyline)	Human Plasma	0.5-200 ng/mL 0.5-200 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Madison	Active	P558.01
	Rabbit Plasma	1-500 ng/mL 1-500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P647.00
Amitriptyline	Human Plasma	0.5-200 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS301
Amlodipine	Human Plasma	0.1-50 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMS151
Amoxicillin	Human Serum	50-25,000 ng/ml	LC/MS/MS	0.1 mL	None	-70	Richmond	Active	LCMS141

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Amoxicillin	Human Plasma	50-25,000 ng/ml	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMS141.2
Atenolol	Human Whole Blood	10-1,000 ng/mL	HPLC/FL	1.0 mL	Sodium Heparin	-20	Madison	Inactive	P161.01
Aripiprazole (Dehydroaripiprazole)	Human Plasma	0.5-250 ng/mL 0.1-50 ng/mL	LC/MS/MS	0.25 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS 386
Atovaquone	Human Plasma	0.025-10 µg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS360
Atorvastatin (Ortho-hydroxy atorvastatin) (Para-hydroxy atorvastatin)	Human Plasma	0.1-75 ng/mL 0.1-75 ng/mL 0.1-75 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Richmond	Active	LCMS120
	Human Plasma	0.1-75 ng/mL 0.1-75 ng/mL 0.1-75 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P913.00
	Rat Plasma	0.1-75 ng/mL 0.1-75 ng/mL 0.1-75 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Richmond	Active	LCMS120.3
Avastin aka Bevacizumab	Human Serum	10-2,560 ng/mL	MSD-ECL	0.025 mL	None	-70	Richmond	Active	ICD343
Azithromycin	Human Plasma	2.5-1,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Inactive	LCMS168
Azidothymidine (AZT)	Human Plasma	50-25,000 ng/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Richmond	Inactive	LC96
	Human Urine	0.2-25 µg/mL	HPLC/UV	1.0 mL	None	-20	Richmond	Inactive	LC96.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Beclomethasone (Beclomethasone-17-propionate) (Beclomethasone dipropionate)	Human Plasma	1-250 pg/mL 1-250 pg/mL 1-250 pg/mL	LC/MS/MS	1.0 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC333
Benazepril (Benazeprilat)	Human Plasma	10-1,000 ng/mL 10-1,000 ng/mL	LC/MS/MS	0.05 mL	Na <sub>2</sub> EDTA	-20	Richmond	Active	LCMS107
Bicalutamide	Human Plasma	10-1,250 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS318
<b>Bisoprolol</b>	<b>Human Plasma</b>	<b>0.5-200 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.1 mL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Richmond</b>	<b>Active</b>	<b>LCMSB501</b>
<b>Bortezomib</b>	<b>Human Plasma</b>	<b>0.5-500 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.1 mL</b>	<b>Sodium Heparin</b>	<b>-70</b>	<b>Richmond</b>	<b>Active</b>	<b>LCMS509</b>
Brompheniramine	Human Plasma	0.05-20 ng/mL	LC/MS/MS	0.1 mL	Lithium Heparin	-20	Richmond	Active	LCMSC436
Budesonide	Human Plasma	10-2,500 pg/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P429.00
Buprenorphine (Norbuprenorphine) & Naloxone	Human Plasma	50-10,000 pg/mL 50-10,000 pg/mL 5-500 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P942.01
Bupropion (Hydroxybupropion) (Erythrohydroxybupropion) (Threoxyhydroxybupropion)	Human Plasma	0.5-200 ng/mL 1-400 ng/mL 1-400 ng/mL 1-400 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-70	Madison	Active	P986.00
Butalbital	Human Plasma	0.05-25 µg/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Richmond	Active	LC196
	Human Plasma	0.1-10 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	<b>Inactive</b>	P301.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Caffeine (1,7-Dimethylxanthine)	Human Plasma	20-20,000 ng/mL 20-20,000 ng/mL	LC/MS/MS	0.15 mL	Sodium Heparin	-20	Richmond	Active	LCMSC207
	Human Urine	1-50 µg/mL 1-50 µg/mL	HPLC/UV	0.1 mL	None	-20	Madison	Active	P590.02
	Human Plasma	25-25,000 ng/mL 25-25,000 ng/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Madison	Active	P589.02
	Human Plasma	25-25,000 ng/mL 25-25,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P734.03
	Human Urine	50-50,000 ng/mL 50-50,000 ng/mL	LC/MS/MS	0.025 mL	None	-20	Madison	Active	P739.03
Candesartan	Human Plasma	0.5-400 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC451
	Human Urine	10-5,000 ng/mL	LC/MS/MS	0.05 mL	None	-20	Richmond	Active	LCMSC451.1
Carbidopa & Levodopa	Human Plasma	0.5-250 ng/mL 5-2,500 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P838.00
Carbidopa (3-O-MethylDopa) & Levodopa	Human Plasma	0.5-250 ng/mL 25-25,000 ng/mL 10-5,000 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P959.00
Carbamazepine (Carbamazepine 10,11 Epoxide)	Human Plasma	0.05-20 µg/mL 0.05-20 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Active	P239.03

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Carbamazepine (Carbamazepine 10,11 Epoxide)	Human Plasma	0.05-50 µg/mL 0.05-50 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P626.00
Carbamazepine, Free (unbound) (Carbamazepine 10,11 Epoxide)	Human Plasma Ultrafiltrate	0.1-20 µg/mL 0.1-20 µg/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Madison	Inactive	P238.02
(Carbamazepine 10,11 Epoxide)	Human Urine	0.1-20 µg/mL	HPLC/UV	0.25 mL	None	-20	Madison	Inactive	P257.01
(Monohydroxycarbamazepine)	Human Plasma	1-50 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-70	Madison	Active	P886.00
Carvedilol	Human Plasma	0.1-50 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P789.01
Cefaclor	Human Plasma	0.1-25 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-70	Richmond	Active	LC305
Cefadroxil	Human Serum	0.25-50 µg/mL	HPLC/UV	0.10 mL	None	-70	Richmond	Inactive	LC52
	Human Urine	5-1,000 µg/mL	HPLC/UV	0.05 mL	None	-20	Richmond	Inactive	LC52.1
Celecoxib	Human Plasma	5-2,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSD417
Celiprolol	Human Plasma	10-2,000 ng/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Richmond	Inactive	LC30
Cephalexin	Human Plasma	0.1-100 µg/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB382
Cetirizine	Human Plasma	0.5-500 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB169
Chlorpheniramine	Human Plasma	0.25-50 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMS214
Chlorthalidone	Human Plasma	2-1,000 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS389.2

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Chlorthalidone	Human Urine	10-5,000 ng/mL	LC/MS/MS	0.1 mL	None	-20	Richmond	Active	LCMS389.1
Chlorzoxazone (6-Hydroxychlorzoxazone)	Human Plasma	5-5,000 ng/mL 5-5,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMS209
Cimetidine	Human Urine	5.0-5,000 µg/mL	HPLC/UV	0.05 mL	None	-20	Richmond	Active	LC358.1
	Human Plasma	0.01-10 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Active	LC358
Ciprofloxacin	Human Plasma	25-10,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS176
Citalopram	Human Plasma	1-500 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMS332
Citalopram (Desmethylcitalopram)	Human Plasma	1-500 ng/mL 0.2-100 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P787.00
Clarithromycin	Human Plasma	20-10,000 ng/mL 5-2,500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB179
Clavulanic Acid	Feline Plasma	0.1-10 µg/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-80	Madison	Active	P997.00
	Canine Plasma	0.1-10 µg/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-70	Madison	Active	P990.00
Clopidogrel (Clopidogrel Acid)	Human Plasma	0.2-200 ng/mL 10-10,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS271
Clorsulon	Bovine Plasma	0.1-10 ng/mL	LC/MS/MS	1.0 mL	Sodium Heparin	-20	Madison	Active	P717.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Cocktail Assay, "Epilepsy" Carbamazepine Lamotrigine Phenobarbital Phenytoin Topiramate	Human Plasma	0.1-50 µg/mL 0.1-20 µg/mL 0.1-50 µg/mL 0.1-30 µg/mL 0.1-50 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P837.00
Cocktail Assay, "Antiviral" Darunavir Lopinavir Raltegravir Ritonavir	Human Plasma	10-10,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P1002.00
Codeine	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Madison	Active	P753.01
Codeine & Morphine (Morphine-3-β-glucuronide) (Morphine-6-β-glucuronide)	Human Plasma	1-100 ng/mL 0.2-20 ng/mL 2-200 ng/mL 0.5-50 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSC396
Codeine & Morphine & Guafenesin	Human Plasma	0.1-100 ng/mL 0.2-20 ng/mL 2-2,000 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Madison	Active	P919.02
Cortisol (6β-Hydroxycortisol)	Human Urine	1-200 ng/mL 6-1,200 ng/mL	LC/MS/MS	0.1 mL	None	-20	Richmond	Active	LCMSC355
Cortisol <i>aka hydrocortisone</i> Cortisone	Human Urine	1-250 ng/mL 1-250 ng/mL	LC/MS/MS	0.05 mL	None Requires CHAPS	-70	Richmond	Active	LCMSC472

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Cortisol <i>aka hydrocortisone</i> Cortisone	Human Serum	1-500 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.05 mL	None	-70	Richmond	Active	LCMSD472.1
5- $\alpha$ -Tetrahydrocortisol (5- $\alpha$ -THF) 5- $\beta$ -Tetrahydrocortisol (5- $\beta$ -THF) Tetrahydrocortisone (THE)	Human Urine	20-5,000 ng/mL 20-5,000 ng/mL 20-5,000 ng/mL	LC/MS/MS	0.05 mL	None <b>Requires CHAPS</b>	-70	Richmond	Active	LCMSC471
Cyclophosphamide	Human Plasma	0.25–50 $\mu$ g/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Richmond	Active	LC232
D4T	Human Plasma	0.025-25 $\mu$ g/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	<b>Inactive</b>	LC211
Danazol	Human Plasma	2-200 ng/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Richmond	<b>Inactive</b>	LC18
<b>Darunavir</b> <b>See "Cocktail Assay Antiviral"</b>	<b>Human Plasma</b>	<b>10-10,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.1 mL</b>	<b>K<sub>3</sub>EDTA</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P1002.00</b>
Daunorubicin (Daunorubicinol)	Human Plasma	5-200 ng/mL 5-200 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Richmond	<b>Inactive</b>	LC301
Desipramine (2-Hydroxydesipramine)	Human Plasma	0.25-100 ng/mL 0.25-100 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMS103
	Human Plasma	0.25-100 ng/mL 0.25-100 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P644.01
Dexamethasone	Human Plasma	0.5-500 ng/mL	LC/MS/MS	0.2 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS328
Dextromethorphan	Human Plasma	10 -5,000 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P732.01
Dextromethorphan (Dextrophan)	Dog Plasma	2.5-500 ng/mL 2.5-500 ng/mL	HPLC/FL	0.2 mL	Sodium Heparin	-20	Madison	Active	P423.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Dextromethorphan (Dextrorphan)	Rat Plasma	5-1,000 ng/mL 5-1,000 ng/mL	HPLC/FL	0.1 mL	Sodium Heparin	-20	Madison	Inactive	P422.00
	Human Plasma	0.05- 50 ng/mL 0.8-800 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P668.05
Dextromethorphan, Total** & Unconjugated (Dextrorphan) (3-methoxymorphinan) (3-hydroxymorphinan)	Human Urine	0.001-1.0 µg/mL 0.020-20 µg/mL 0.001-1.00 µg/mL 0.01-10.0 µg/mL	LC/MS/MS	0.05 mL	None	-20	Madison	Active	P600.07
Dextromethorphan & Pseudoephedrine	Human Plasma	10-5,000 pg/mL 0.8-400 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P617.02
Diazepam (N-Desmethyldiazepam)	Human Plasma	2-1,000 ng/mL 1-500 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS 466
Diclofenac	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSC368
Diclofenac (4'-hydroxydiclofenac)	Human Plasma	10-2,000 ng/mL 2.5-500 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMSB368.1
Digoxin	Human Plasma	0.15-8 ng/mL	RIA	0.2 mL	Sodium Heparin	-70	Richmond	Active	ICD11.1
	Human Serum	0.15-8 ng/mL	RIA	0.1 mL	None	-70	Richmond	Active	ICD11.3
	Human Urine	1-40 ng/mL	RIA	0.1 mL	None	-70	Richmond	Active	ICD11.2
	Human Plasma	0.1-4 ng/mL	RIA	0.4 mL	EDTA	-70	Richmond	Active	ICD11.5
Digoxin	Human Plasma	0.1-50 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC455.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Diltiazem (Desacetyldiltiazem) (Desmethyldiltiazem)	Human Plasma	0.1-200 ng/mL 0.1-200 ng/mL 0.1-200 ng/mL	LC/MS/MS	0.4 mL	Sodium Heparin	-20	Madison	Active	P672.03
Docetaxel	Human Plasma	10-5,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB216
Donepezil	Human Plasma	0.2-50 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMSB203
	Human Plasma	0.2-100 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P983.00
Doxazosin	Human Plasma	0.05-50 ng/mL	HPLC/FL	1.0 mL	Sodium Heparin	-20	Richmond	<b>Inactive</b>	LC266
Doxepin (Desmethyldoxepin)	Human Plasma	0.5-250 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P815.01
Doxorubicin	Human Plasma	1-200 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-70	Richmond	<b>Inactive</b>	LC252
Duloxetine	Human Plasma	0.2-200 ng/mL	LC/MS/MS	0.5 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P882.00
Efavirenz	Human Plasma	0.10-10 µg/mL	HPLC/UV	0.1 mL	Sodium Heparin	-20	Richmond	Active	LC346
Emtricitabine & Tenofovir	Human Plasma	20-4,000 ng/mL 5-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC 468
	Human Plasma	5-1,000 ng/mL 20-4,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P947.01
Enalapril (Enalaprilat)	Human Plasma	0.25-100 ng/mL 0.25-100 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Madison	Active	P435.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Enalapril (Enalaprilat)	Human Plasma	0.5-250 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Inactive	LCMS31
	Dog Plasma	0.25-400 ng/mL 0.25-400 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Madison	Active	P532.01
Epinephrine	Human Plasma	0.05-8 ng/mL	HPLC/FL	1.0 mL	Sodium Metabisulfite	-20	Madison	Active	P981.02
Erlotinib	Human Plasma	10-10,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSC390
	<b>Human Plasma</b>	<b>5-5,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.025 mL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P1008.00</b>
Erythromycin	Human Plasma	0.2-20 µg/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS463
Estazolam	Human Plasma	0.1-30 ng/mL	LC/MS/MS	1.0 mL	Sodium Heparin	-20	Madison	Inactive	P336.00
β-Estradiol & Equilin & Estrone	Human Plasma	5.0-500 pg/mL 5.0-500 pg/mL 5.0-500 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-70	Richmond	Active	LCMSC233
	Human Plasma	1.0-100 pg/mL 5.0-500 pg/mL 5.0-500 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC350
	Human Serum	1.0-100 pg/mL 5.0-500 pg/mL 5.0-500 pg/mL	LC/MS/MS	0.5 mL	None	-70	Richmond	Active	LCMSC350.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
β-Estradiol (Estrone)	Human Serum	2.5 – 250 pg/mL 5.0 – 500 pg/mL	LC/MS/MS	0.5 mL	None	-20	Richmond	Inactive	LCMSC248.1
	Human Plasma	2.5 – 250 pg/mL 5 – 500 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC248
β-Estradiol, Total**	Human Serum	0.025-10 ng/mL	LC/MS/MS	0.1 mL	None	-70	Richmond	Active	LCMS295
Estrone Sulfate	Human Serum	50-5,000 pg/mL	LC/MS/MS	0.5 mL	None	-70	Richmond	Active	LCMS27.1
Estrone & Equilin, Total**	Human Plasma	10-5,000 pg/mL 10-5,000 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC351
	Human Serum	10-5,000 pg/mL 10-5,000 pg/mL	LC/MS/MS	0.5 mL	None	-70	Richmond	Active	LCMSC351.1
Ethanol	Human Plasma	20-1,000 µg/mL	GC/FID	0.1 mL	Potassium Oxalate/Sodium Fluoride	-20	Richmond	Active	GC89
Ethinyl Estradiol & Norethindrone	Human Plasma	2-500 pg/mL 50-25,000 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Richmond	Active	LCMSC255
Ethinyl Estradiol & Norgestrel	Human Plasma	2-500 pg/mL 50-25,000 pg/mL	LC/MS/MS	0.5 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC256
Ethinyl Estradiol Sulfate	Human Plasma	50-10,000 pg/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSC398
Etodolac	Human Plasma	0.2-50 µg/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Richmond	Active	LC193
Etoricoxib	Human Plasma	5-2500 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P971.01

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Ezetimibe	Human Plasma	0.2-100 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC439
	Human Plasma	0.2-100 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P945.00
Ezetimibe, Total	Human Plasma	1.0-500 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC440
	Human Plasma	1.0-500 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P944.00
Famotidine	Human Plasma	5-250 ng/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Madison	Inactive	P255.01
Fenofibric Acid	Human Plasma	0.05-30 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSB235
Fentanyl  Fentanyl (Norfentanyl)	Human Plasma	25-10,000 pg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSB186.5
	Human Plasma	5-2,000 pg/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC186.4
	Human Serum	10-2,000 ng/mL	LC/MS/MS	0.1 mL	None	-20	Madison	Active	P869.01
	Dog Plasma	40-40,000 pg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P966.00
	Mouse Plasma	0.04-40 ng/mL 0.04-40 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P637.02
	Rat Plasma	0.04-40 ng/mL 0.04-40 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P823.01
	Rabbit Plasma	0.04-40 ng/mL 0.04-40 ng/mL	LC/MS/MS	0.2 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P636.01
Fentanyl (Norfentanyl) (DespropionylFentanyl)	Human Plasma	10-500 pg/mL 10-500 pg/mL 10-500 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P624.02

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Fexofenadine	Human Plasma	0.5-500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P745.04
	Human Urine	0.05-10 µg/mL	LC/MS/MS	0.02 mL	None	-20	Madison	Active	P746.04
Fluconazole	Human Plasma	20-10,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P791.00
Flunixin	Equine Plasma	10-5,000 ng/mL	LC/MS/MS	0.025 mL	Sodium Heparin	-70	Madison	Active	P970.00
Fluoxetine (Norfluoxetine)	Human Plasma	0.25-100 ng/mL 0.25-100 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMS52
Fluticasone Propionate	Human Plasma	2.5-500 pg/mL	LC/MS/MS	1.0 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC173
	Human Plasma	1.0-500 pg/mL	LC/MS/MS	1.0 mL	Potassium Oxalate/Sodium Fluoride	-20	Madison	Active	P817.02
Fluvoxamine	Human Plasma	0.5-500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P658.01
Formoterol	Human Plasma	0.5-250 pg/mL	LC/MS/MS	1.0 mL	Lithium Heparin with Eserine hemisulfate	-20	Madison	Active	P860.02
	Human Plasma	0.25-125 pg/mL	LC/MS/MS	1.0 mL	Lithium Heparin with Eserine hemisulfate	-20	Madison	Active	P996.00
	Human Urine	6.0-1,200 pg/mL	LC/MS/MS	0.1 mL	None	-20	Madison	Active	P883.01
Formoterol Enantiomers (R,R and S,S)	Human Urine	6.0-1,200 pg/mL each	LC/MS/MS	0.5 mL	None	-20	Madison	Active	P868.00
Furosemide	Human Urine	0.05-15 µg/mL	HPLC/FL	1.0 mL	None	-20	Madison	Active	P376.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Human Plasma	0.02-10 µg/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P372.01
	Human Plasma	0.02-10 µg/mL	HPLC/FL	0.2 mL	Sodium Heparin	-20	Richmond	Active	LC354
	Human Urine	0.2-50 µg/mL	LC/MS/MS	0.2 mL	None	-20	Richmond	Active	LCMS317
Gabapentin	Human Plasma	0.05- 10 µg/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS359
Galanthamine	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB312
Gatifloxacin	Human Plasma	0.025-10 µg/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB285
Gemfibrozil	Human Plasma	0.05-25 µg/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC56
	Human Plasma	0.05-50 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin and EDTA	-20	Madison	Inactive	P115.00
Glimepiride	Human Plasma	2-1,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS266
Glipizide	Human Plasma	1-500 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS137
Glyburide	Human Plasma	1-500 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS159
Glyburide <i>Lower range for pediatric studies</i>	Human Plasma	0.10-50 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS159.1
<b>Guaifenesin</b>	<b>Human Plasma</b>	<b>2-2,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.075 mL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P973.00</b>
Guaifenesin & Codeine & Morphine	Human Plasma	2-2,000 ng/mL 0.1-100 ng/mL 0.2-20 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Madison	Active	P919.02

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Herceptin <i>aka Trastuzumab</i>	Human Serum	10-200 µg/mL	ELISA	0.2 mL	None	-70	Richmond	Active	ICD178
	Human Plasma	10-100 µg/mL	ELISA	0.1 mL	EDTA	-20	Richmond	Active	ICD178.1
Hydrochlorothiazide	Human Plasma	2-200 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSB205.4
	Human Plasma	2-200 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB205
	Human Urine	0.1-25 µg/mL	LC/MS/MS	0.05 mL	None	-20	Richmond	Active	LCMSB205.3
Hydrocodone (Hydromorphone)	Human Plasma	0.1-100 ng/mL 0.05-50 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC492
	Human Plasma	0.1-100 ng/mL 0.05-50 ng/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Inactive	P623.00
	Hamster Plasma	1-400 ng/mL 1-400 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Inactive	P659.00
Hydrocodone & Ibuprofen	Human Plasma	0.25-50 ng/mL 0.25-50 µg/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P774.01
Ibuprofen	Human Plasma	0.2-50 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB409
	Human Plasma	0.2-50 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Inactive	P046.04
IFN Gamma	Human Serum	31.3-1,000 pg/mL	ELISA	0.2 mL	None	-70	Richmond	Active	ICD278
IL-1β	Human Serum	7.82-250 pg/mL	ELISA	0.4 mL	None	-20	Richmond	Active	ICD151

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
IL-2	Human Serum	62.5-2,000 pg/mL	ELISA	0.2 mL	None	-20	Richmond	Active	ICD150
IL-4	Human Serum	31.3-2,000 pg/mL	ELISA	0.1 mL	None	-70	Richmond	Active	ICD277
IL-6	Human Serum	5-320 pg/mL	ELISA	0.2 mL	None	-20	Richmond	Active	ICD141
IL-10	Human Serum	10- 500 pg/mL	ELISA	0.3 mL	None	-20	Richmond	Active	ICD142
Irbesartan	Human Plasma	2.5-2,500 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB232
Isopropamide	Human Plasma	0.05-10 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Inactive	P409.01
Itraconazole (2-Hydroxyitraconazole)	Human Plasma	0.5-100 ng/mL 0.5-100 ng/mL	HPLC/FL	1.0 mL	Sodium Heparin	-20	Madison	Active	P476.01
Ivermectin B1a	Equine Plasma	0.5–50 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P575.01
	Equine Serum	0.50-50 ng/mL 0.84-1.68 ng/mL	HPLC/FL	2.0 mL	None	-20	Madison	Active	P389.00
Ivermectin B1b	Equine Plasma	0.10-10 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P621.00
Ivermectin B1a	Swine Plasma	0.50-100 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P646.00
	Bovine Plasma	0.50-100 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P479.01
	Dog Plasma	0.50-50 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P457.02
Ivermectin B1a	Dog Plasma	0.10-10 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Madison	Active	P715.01
Ketamine (Norketamine)	Human Serum	5-5,000 ng/mL 2.5-2,500 ng/mL	LC/MS/MS	0.2 mL	None	-20	Madison	Active	P627.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Rabbit Plasma	10-10,000 ng/mL 5-5,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P648.00
	Human Plasma	0.5-500 ng/mL 0.5-500 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P663.02
<b>S-Ketamine and R-Ketamine S-Norketamine and R-Norketamine</b>	<b>Feline Plasma</b>	<b>20-8,000 ng/mL 10-4,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>50 uL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P703.01</b>
Ketoconazole	Human Plasma	0.05-5 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB202
Ketoprofen	Human Plasma	1-500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB162
	Human Plasma	0.1-50 ng/mL	LC/MS/MS	0.4 mL	Sodium Heparin	-20	Madison	Active	P747.02
	Mini Pig	0.1-50 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P821.00
Ketoprofen, and Total**	Human Plasma	1.0-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P767.00
	Human Urine	5-5,000 ng/mL	LC/MS/MS	0.025 mL	None	-20	Madison	Active	P768.01
Ketoprofen Enantiomers	Human Plasma	0.5– 250 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P782.00
Ketorolac	Human Plasma	0.025-25 µg/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Richmond	Active	LC207
Lamivudine (3TC)	Human Plasma	10-4000 ng/mL	HPLC/UV	0.5 mL	EDTA	-20	Madison	Active	P560.00
	Human Plasma	10-5,000 ng/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Active	LC269
Lamotrigine See Cocktail Assay "Epilepsy"	Human Plasma	0.1-20 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P837.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Lansoprazole	Human Plasma	2-2,000 ng/mL	LC/MS/MS	0.05mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS178
	Human Plasma	2-2,000 ng/mL	LC/MS/MS	0.05mL	K <sub>3</sub> EDTA	-20	Madison	Active	P896.02
Leucovorin Diastereomers (6S-5-Methyltetrahydrofolic Acid Diastereomers)	Human Serum	0.5-250 µg/mL 0.1-50 µg/mL	LC/MS/MS	0.05 mL	None	-70	Richmond	Active	LCMSC387
Leuprolide	Human Plasma	0.025-25 ng/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC245
Levetiracetam	Human Plasma	0.05-50 µg/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P811.03
Levodopa & Carbidopa	Human Plasma	5-2,500 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P838.00
Levodopa & Carbidopa (3-O-MethylDopa)	Human Plasma	10-5,000 ng/mL 0.5-250 ng/mL 25-25,000 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P959.00
Levofloxacin	Human Serum	25-15,000 ng/mL	HPLC/UV	0.5 mL	None	-20	Richmond	Active	LC287
Lidocaine	Human Plasma	2-500 ng/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Madison	Active	P654.00
	Pig Plasma	2-500 ng/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Madison	Active	P634.00
	Rat Plasma	1-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P667.00
Lidocaine	Human Plasma	0.2-200 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P871.00
Lisinopril	Human Plasma	0.5-400 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Madison	Active	P698.01
	Human Urine	25-25,000 ng/mL	LC/MS/MS	0.025 mL	None	-20	Madison	Active	P694.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Rat Plasma	1.0-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P673.00
Loperamide	Human Plasma	10-10,000 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P557.02
	Human Plasma	10-10,000 pg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSC310
<b>Lopinavir See Cocktail Assay "Antiviral"</b>	<b>Human Plasma</b>	<b>10-10,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.1 mL</b>	<b>K<sub>3</sub>EDTA</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P1002.00</b>
Loratadine (Descarboethoxyloratadine) (3-OHDescarboethoxyloratadine)	Human Plasma	0.025-10 ng/ml 0.025-10 ng/ml 0.025-10 ng/ml	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS188.1
Lorazepam	Human Plasma	0.5-50 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMS142
Lorazepam, Total**	Human Plasma	2.5-500 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSB294
Losartan (Acid Metabolite)	Human Plasma	0.5-500 ng/mL 0.5-500 ng/mL	LC/MS/MS	0.15 mL	Sodium Heparin	-20	Richmond	Active	LCMSC208
	Human Urine	5-2,000 ng/mL 2-2,000 ng/mL	LC/MS/MS	0.1 mL	None	-20	Richmond	Inactive	LCMSC208.1
	Human Plasma	0.05-12.5 ng/mL 0.1-25 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-70	Richmond	Inactive	LCMS267
Meclocycline	Human Plasma	25-10,000 ng/mL	HPLC/UV	0.2 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LC337
Medroxyprogesterone Acetate	Human Plasma	0.01-10 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LCMS99
	Human Serum	0.02-5 ng/mL	LC/MS/MS	0.25 mL	None	-20	Richmond	Active	LCMSC198
Memantine	Human Plasma	0.5-50 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P786.01

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Human Plasma	0.5-50 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P786.01
Metformin	Human Plasma	2-2,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMS153
	Human Plasma	2-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS153.5
	Human Urine	10-10,000 µg/mL	HPLC/UV	0.05 mL	None	-20	Richmond	Active	LC251.1
	Dog Plasma	1-1,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB153.4
Methadone	Human Plasma	0.1-50 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS323
	Human Urine	25-1,000 ng/mL	HPLC/UV	1.0 mL	None	-20	Richmond	Inactive	LC238.2
Methadone Enantiomers	Human Plasma	5-1,000 ng/mL 5-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC323.1
Methotrexate (7-Hydroxymethotrexate)	Human Plasma	5-500 ng/mL 5-500 ng/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Madison	Active	P461.02
	Human Urine	0.1-10 ug/mL 0.1-10 ug/mL	HPLC/UV	0.05 mL	None	-20	Madison	Active	P683.02
Methotrexate (7-Hydroxymethotrexate)	Human Plasma	1.0-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSC380.1
	Human Urine	0.05-25 ug/mL	LC/MS/MS	0.05 mL	None	-20	Richmond	Active	LCMSC380.2
	Human Plasma	1.0-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMS380
Methylphenidate	Human Plasma	0.25-25 ng/mL	LC/MS/MS	0.10 mL	Sodium Heparin	-20	Richmond	Inactive	LCMSB187

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Methylphenidate (Ritalinic Acid)	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMS230.1
	Beagle Plasma	0.10-100 ng/mL 2.0-2,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS230
	Human Plasma	0.1-50 ng/mL 2.0-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P689.02
Methylphenidate Enantiomers (Ritalinic Acid Enantiomers)	Human Plasma	0.05-50 ng/mL (D) 0.01-10 ng/mL (L) 1.0-1,000 ng/ml (D) 1.0-1,000 ng/ml (L)	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P779.00
Methyltestosterone	Human Plasma	5- 25,000 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC349
Metoprolol	Human Plasma	5-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P940.00
Metronidazole (1-Hydroxymetronidazole)	Human Plasma	10-1,000 ng/mL 10-1,000 ng/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Active	P686.00
Mevalonic Acid	Human Plasma	0.5-20 ng/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P758.02
	Human Urine	25-1,000 ng/mL	LC/MS/MS	0.25 mL	None	-20	Madison	Active	P759.00
Midazolam (1-Hydroxy midazolam) (1-HydroxyMidazolam, Total**)	Human Plasma	0.1-100 ng/mL 0.1-100 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Richmond	Active	LCMS97
	Human Urine	1-1,000 ng/mL	LC/MS/MS	0.2 mL	None	-20	Richmond	Inactive	LCMS228
Midazolam (1-Hydroxy midazolam) (4-Hydroxy midazolam)	Human Plasma	0.1-100 ng/mL 0.1-100 ng/mL 0.1-100 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P594.05

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Midazolam (1-HydroxyMidazolam)	Human Urine	0.05-50 ng/mL 1-1,000 ng/mL	LC/MS/MS	0.2 mL	None	-20	Madison	Active	P678.04
	Human Urine	0.05-50 ng/mL 2-2,000 ng/mL	LC/MS/MS	0.2 mL	None	-20	Madison	Active	P870.01
Minocycline	Human Plasma	0.02-10 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Active	LC63
	Human Serum	0.02-10 µg/mL	HPLC/UV	0.5 mL	None	-20	Richmond	Active	LC63.1
	Human Saliva	0.02-10 µg/mL	HPLC/UV	0.5 mL	None	-20	Richmond	Active	LC63.5
	Human Crevicular Fluid	10-2,000 ng/µL	HPLC/UV	1.0 µL	None	-20	Richmond	Active	LC360
Mirtazapine	Human Plasma	0.05-100 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMS145
Mometasone Furoate	Human Plasma	1-250 pg/mL	LC/MS/MS	1.0 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSD345
	Human Plasma	0.25-25 pg/mL	LC/MS/MS	1.0 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSD345.1
	Dog Plasma	10-1,000 pg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSD345.4
Mometasone Furoate	Rat Plasma	10-1,000 pg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSD345.3
Montelukast	Human Plasma	5-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-80	Richmond	Active	LCMSC385
Morphine (Morphine-3-β-glucuronide) (Morphine-6-β-glucuronide)	Human Plasma	0.5-50 ng/mL 10-1,000 ng/mL 2-200 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS 454

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Morphine (Morphine-3-β-glucuronide) (Morphine-6-β-glucuronide) Codeine	Human Plasma	0.2-20 ng/mL 2-2,00 ng/mL 0.5-50 ng/mL 1-100 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSC396
Morphine & Codeine & Guafenesin	Human Plasma	0.2-20 ng/mL 0.1-100 ng/mL 2-2,000 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Madison	Active	P919.02
Moxifloxacin	Human Plasma	25-5,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSB276
Naloxone & Buprenorphine (Norbuprenorphine)	Human Plasma	5-500 pg/mL 50-10,000 pg/mL 50-10,000 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P942.01
Naloxone, Total	Human Plasma	25-10,000 pg/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P943.00
Naproxen	Human Plasma	0.1-100 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSB486
Navelbine	Human Plasma	2-500 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-70	Richmond	Inactive	LC253.1
Nevirapine	Human Plasma	25-10,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS250
Niacin <i>aka Vitamin B3, or nicotinic acid</i> (Nicotinuric acid) <i>aka NUA</i>  (1-Methyl-Nicotinamide) <i>aka NMA</i> (N-Methyl-2-Pyridone-5-carboxamide) <i>aka 2PY</i>	Human Plasma	5-2,000 ng/mL 5-2,000 ng/mL	LC/MS/MS	0.25mL	Sodium Heparin	-70	Richmond	Active	LCMSB264
	Human Urine	0.02 -20 µg/mL 0.2 - 200 µg/mL	LC/MS/MS	0.1 mL	None	-20	Richmond	Active	LCMSB264.1
	Human Urine	0.5-250 µg/mL 2.5-250 µg/mL	LC/MS/MS	0.05 mL	None	-20	Richmond	Active	LCMS274.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Human Plasma	0.1-10 µg/mL 0.2- 20 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB274
	Human Plasma	0.1-10 µg/mL 0.2- 20 µg/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSB274.2
Nicotine (Cotinine)	Human Plasma	0.5-100 ng/mL 2.5-500 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P538.02
Nizatidine (N-Desmethylnizatidine)	Human Plasma	5-5,000 ng/mL 5-5,000 ng/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Richmond	Active	LC264
Norelgestromin <i>aka 17-desacetyl norgestimate</i>	Human Plasma	0.02-10 ng/mL	LC/MS/MS	0.5 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC372
Norethindrone & Ethinyl Estradiol	Human Plasma	50-25,000 pg/mL 2-500 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Richmond	Active	LCMSC255
Norgestrel & Ethinyl Estradiol	Human Plasma	50-25,000 pg/mL 2-500 pg/mL	LC/MS/MS	0.5 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC256
Octreotide	Human Plasma	0.025-25 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC296
Octreotide	Monkey Plasma	5.0-250 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC296.2
	Monkey Plasma	0.1-25 ng/mL	LC/MS/MS	0.2 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC296.1
	Rabbit Plasma	1-100 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC296.3
Olanzapine (Desmethyloanzapine)	Human Plasma	0.05-50 ng/mL 0.1-100 ng/mL	LC/MS/MS	0.3 mL	Sodium Heparin	-20	Richmond	Active	LCMS84

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Omeprazole	Human Plasma	1-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB131
Omeprazole (5-Hydroxyomeprazole) (Omeprazole sulphone)	Human Plasma	1-1,000 ng/mL 1-1,000 ng/mL 1-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P749.04
Omeprazole (5-Hydroxyomeprazole)	Equine Plasma	1-1,000 ng/mL 1-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P822.00
Ondansetron	Human Plasma	0.5-250 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS165
	Dog Plasma	0.5-250 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS165.1
Oxaprozin	Human Plasma	0.5-500 µg/mL	HPLC/UV	0.1 mL	Sodium Heparin	-20	Richmond	Inactive	LC243
Oxcarbazepine (10-Hydroxycarbamazepine (MHD)) (10,11-Dihydroxycarbamazepine (DHD))	Human Plasma	0.1-5 µg/mL 0.1-5 µg/mL 0.5-25 µg/mL	HPLC/UV	0.3 mL	Sodium Heparin	-70	Madison	Active	P743.00
Oxcarbazepine	Human Plasma	0.1-5 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P850.00
(Monohydroxyoxcarbazepine <i>aka</i> MHD)	Human Plasma	1-50 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-70	Madison	Active	P886.00
Oxybutynin Enantiomers (Desethyloxybutynin enantiomers)	Human Plasma	0.5-200 ng/mL 0.5-200 ng/mL	LC/MS/MS	0.2 mL	K <sub>2</sub> EDTA	-70	Richmond	Active	LCMS 411
Oxycodone	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P649.01
Oxymorphone	Human Plasma	10-1,000 pg/mL	LC/MS/MS	0.5 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P820.01

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Paclitaxel (Taxol ®)	Human Plasma	5-1,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Richmond	Active	LCMS163
	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS163.2
Paroxetine	Human Plasma	0.25-50 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Richmond	Active	LCMS161
	Human Plasma	0.25-50 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin K <sub>3</sub> EDTA	-20	Madison	Active	P776.01
Pemoline	Human Plasma	25-5,000 ng/mL	HPLC/UV	0.2 mL	Sodium Heparin	-20	Richmond	Inactive	LC258
Pentoxyfilline (Hydroxypentifylline) (Carboxypropylpentifylline)	Human Plasma	5-1,000 ng/mL 10-2,000 ng/mL 10-2,000 ng/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Madison	Active	P348.01
<b>Phenethylamine</b>	<b>Human Urine</b>	<b>1.13-501 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.2 mL</b>	<b>None</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P655.00</b>
Phenobarbital See Cocktail Assay "Epilepsy"	Human Plasma	0.1-50 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P837.00
Phenylbutazone	Equine Plasma	0.1-40 µg/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Madison	Active	P554.04
Phenylephrine, Total**	Human Plasma	5.0-1,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-70	Madison	Active	P806.01
	Human Plasma	1.0-1,000 ng/mL	LC/MS/MS	0.05 mL	Lithium Heparin	-70	Richmond	Active	LCMS257
	Human Urine	100-5,000 ng/mL	LC/MS/MS	0.05 mL	None	-70	Madison	Active	P958.01
Phenylephrine (Phenylephrine Glucuronide)  (3-Hydroxymandelic acid)	Human Urine	2.0-200 ng/mL 2.0-200 ng/mL	LC/MS/MS	0.2 mL	None	-70	Madison	Active	P956.00
	Human Urine	100-5,000 ng/mL	LC/MS/MS	0.05 mL	None	-70	Madison	Active	P957.00
Phenylephrine	Human Plasma	10-2,500 pg/mL	LC/MS/MS	0.25 mL	Lithium Heparin	-70	Richmond	Active	LCMSC392.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
	Human Plasma	10-2,500 pg/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P898.01
Phenylpropanolamine	Human Plasma	5-1,000 ng/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC181
	Human Plasma	5-500 ng/mL	HPLC/UV	1.0 mL	EDTA	-20	Madison	Inactive	P002.05
Phenytoin	Human Plasma	20-5,000 ng/mL	LC/MS/MS	0.25 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSB418
Phenytoin Phenytoin, Total **	Human Plasma	0.05-15 ug/mL 0.1- 30 ug/mL	LC/MS/MS	0.05 mL 0.025 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P639.00
Phenytoin, Free (unbound) (5-Hydroxyphenytoin)	Human Plasma Ultrafiltrate	0.05-15 µg/mL 0.05-15 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P742.03
Phenytoin, Total** (5-Hydroxyphenytoin)	Human Plasma	0.05-15 µg/mL 0.1- 30 µg/mL	LC/MS/MS	0.025 mL	Sodium Heparin	-20	Madison	Active	P742.03
	Human Plasma	0.1-20 µg/mL	HPLC/UV	1.0 mL	Sodium Heparin	-20	Madison	Active	P281.00
Phenytoin (5-Hydroxyphenytoin)	Human Urine	0.1-40 µg/mL 0.1-40 µg/mL	HPLC/UV	0.5 mL	None	-20	Madison	Active	P642.01
	(5-Hydroxyphenytoin), Total**	Human Urine	1-200 µg/mL	HPLC/UV	0.05 mL	None	-20	Madison	Active
Pravastatin	Human Plasma	0.5-250 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS47
Praziquantel	Dog Plasma	10-4,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P989.00
	Feline Plasma	10-4,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P831.01

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Prednisone (Prednisolone)	Human Plasma	0.5-1,000 ng/mL 0.5-1,000 ng/mL	LC/MS/MS	0.3 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P955.00
	Rat Plasma	0.25-250 ng/mL 0.25-250 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-70	Madison	Active	P984.00
Pregabalin	Human Plasma	20-20,000 ng/mL	LC/MS/MS	0.5 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC487
Procainamide (N-Acetylprocainamide)	Human Plasma	0.025-10 µg/mL 0.025-10 µg/mL	HPLC/UV	0.25 mL	Sodium Heparin	-20	Richmond	Inactive	LC306
Progesterone	Human Plasma	0.05-50 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-70	Richmond	Active	LCMSC431
	Rabbit Plasma	0.5-50 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-70	Richmond	Active	LCMSC431.1
Propranolol (4-Hydroxypropranolol)	Human Plasma	0.5-250 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Richmond	Active	LCMSB191
	Human Urine	1-250 ng/mL 1-250 ng/mL	LC/MS/MS	0.25 mL	None	-70	Richmond	Active	LCMSB191.1
Propranolol, Total** (4-Hydroxypropranolol)	Human Plasma	2-600 ng/mL 2-600 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-70	Richmond	Active	LCMSB281.2
Pseudoephedrine	Human Plasma	2.5-500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Richmond	Active	LCMSB201
Pseudoephedrine & Dextromethorphan	Human Plasma	0.8-400 ng/mL 10-5,000 pg/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P617.02
Quetiapine	Human Plasma	1-1000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSB430

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Quinidine	Human Plasma	25-5,000 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC330.1
<b>Raltegravir See Cocktail Assay "Antiviral"</b>	<b>Human Plasma</b>	<b>10-10,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.1 mL</b>	<b>K<sub>3</sub>EDTA</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P1002.00</b>
Ramipril (Ramiprilat)	Human Plasma	0.05-50 ng/mL 0.05-50 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMS311
Ranitidine	Human Plasma	2.5-2,500 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P873.01
Ranolazine	Human Plasma	2-2,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P934.01
Repaglinide	Human Plasma	0.2-200 ng/mL	LC/MS/MS	0.1 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS460
<b>Ribavirin</b>	<b>Human Plasma</b>	<b>5-5,000 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.05 mL</b>	<b>K<sub>2</sub>EDTA</b>	<b>-20</b>	<b>Richmond</b>	<b>Active</b>	<b>LCMS510</b>
Rifabutin (25-O-Desacetylriofabutin)	Human Plasma	2-800 ng/mL 2-800 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P587.06
Rifampin	Human Plasma	0.1-10 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P941.02
Risperidone (9-Hydroxyrisperidone)	Human Urine	1-300 ng/mL 1-300 ng/mL	LC/MS/MS	0.2 mL	None	-20	Madison	Active	P362.01
	Human Plasma	0.1-60 ng/mL 0.1-60 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P528.01
	Human Serum	0.1-30 ng/mL 0.1-30 ng/mL	LC/MS/MS	1.0 mL	None	-20	Madison	Active	P361.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Ritonavir	Human Plasma	10-10,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-70	Richmond	Active	LCMS210.2
See Cocktail Assay "Antiviral"	Human Plasma	10-10,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P1002.00
Rivaroxaban	Human Plasma	0.5-500 ng/mL	LC/MS/MS		K <sub>2</sub> EDTA	-20	Madison	Active	P1014.00
Rizatriptan	Human Plasma	0.2-200 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P845.00
Ropinirole	Human Plasma	10-5,000 pg/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P736.01
Rosiglitazone (N-Desmethylrosiglitazone)	Human Plasma	1-500 ng/mL 1-500 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS322.1
Rosiglitazone	Human Plasma	2-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS322
	Rat Plasma	2-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS322.2
Rosuvastatin	Human Plasma	0.1-100 ng/mL	LC/MS/MS	0.2 mL	K <sub>2</sub> EDTA	-70	Richmond	Active	LCMSC393
Rosuvastatin (N-Desmethylrosuvastatin)	Human Plasma	0.1-100 ng/mL 0.05-50 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-70	Madison	Active	P922.03
(Salicylic Acid)	Human Plasma	1-50 µg/mL	HPLC/UV	0.2 mL	Lithium Heparin	-20	Madison	Active	P383.00
	Dog Plasma	1-50 µg/mL	HPLC/UV	0.2 mL	Sodium Heparin	-70	Richmond	Inactive	LC261.1
Acetylsalicylic Acid (Salicylic Acid)	Human Plasma	0.02-10 µg/mL 0.1-50 µg/mL	LC/MS/MS	0.05 mL	Potassium Oxalate/Sodium Fluoride	-70	Richmond	Active	LCMSC425

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Salmeterol	Human Plasma	1-100 pg/mL	LC/MS/MS	1.0 mL	Potassium Oxalate/Sodium Fluoride	-20	Madison	Active	P798.00
Selegiline (N-Desmethylselegiline) (Amphetamine) (Methamphetamine) (Phenethylamine)	Human Plasma	0.01-5 ng/mL 0.04-20 ng/mL 0.05-25 ng/mL 0.05-25 ng/mL 0.01-5 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P985.00
Sertraline	Human Plasma	0.2-80 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P570.01
Sertraline (Desmethylsertraline)	Human Plasma	0.5-200 ng/mL 0.5-200 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMS98.1
Sildenafil (Desmethylsildenafil)	Human Plasma	1-1,000 ng/mL 1-1,000 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin & K <sub>3</sub> EDTA	-70	Richmond	Active	LCMS231
Simvastatin (Simvastatin Acid)	Human Plasma	0.1-50 ng/mL 0.1-50 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-70	Richmond	Active	LCMSC367
Simvastatin (Simvastatin Acid)	Human Plasma	0.1-50 ng/mL 0.1-50 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P980.00
Solifenacin	Human Plasma	0.1-25 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-20	Madison	Active	P950.00
Sotalol	Human Plasma	10-2,000 ng/mL	HPLC/FL	0.5 mL	EDTA	-20	Madison	Active	P403.01
Sparfloxacin	Human Plasma	0.025-20 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Active	LC332
Sufentanil	Human Plasma	1.0-250 pg/mL	LC/MS/MS	0.5 mL	K <sub>2</sub> EDTA	-70	Richmond	Active	LCMSC239
	Human Urine	5.0-1,000 pg/mL	LC/MS/MS	0.25 mL	None	-70	Richmond	Active	LCMSC239.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Sumatriptan	Human Plasma	0.2-100 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS174
	Pig Plasma	0.2-100 ng/mL	LC/MS/MS	0.25 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS174.3
Sunitinib (N-Desethylsunitinib)	Human Plasma	0.1-50 ng/mL 0.1-50 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC457
T3	Human Serum	0.25-8 ng/mL	RIA	0.2 mL	None	-20	Richmond	Active	ICD86
T4	Human Serum	1-32 µg/dL	RIA	0.05 mL	None	-20	Richmond	Active	ICD87
Tacrolimus	Human Whole Blood	0.25-100 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS356
Tamsulosin	Human Plasma	0.05-25 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P682.02
	<b>Dog Plasma</b>	<b>0.05-25 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.5 mL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P750.01</b>
	<b>Pig Plasma</b>	<b>0.05-25 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.5 mL</b>	<b>Sodium Heparin</b>	<b>-20</b>	<b>Madison</b>	<b>Active</b>	<b>P969.00</b>
Telmisartan	Human Plasma	2-1,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC408
Tenofovir & Emtricitabine	Human Plasma	5-1,000 ng/mL 20-4,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMSC 468
	Human Plasma	5-1,000 ng/mL 20-4,000 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P947.00
Terfenadine (Fexofenadine)	Human Plasma	0.05-50 ng/mL 0.5-500 ng/mL	LC/MS/MS	0.5 mL	Sodium Heparin	-20	Madison	Active	P473.01
Terazosin	Human Plasma	0.5-200 ng/mL	HPLC/FL	0.3 mL	Sodium Heparin	-20	Richmond	<b>Inactive</b>	LC270

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Terconazole	Balb/c Mouse Plasma	1-100 ng/mL	LC/MS/MS	0.025 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMSC217.1
Testosterone (Dihydrotestosterone)	Human Serum	50–10,000 pg/mL 50 -10,000 pg/mL	LC/MS/MS	0.15 mL	None	-20	Richmond	Active	LCMSC260.6
Testosterone, Free (unbound)	Human Serum Ultrafiltrate	1-200 pg/mL	LC/MS/MS	0.2 mL	None	-20	Richmond	Active	LCMSC247
Theophylline	Human Plasma	0.05-25 µg/mL	HPLC/UV	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC9
	Human Serum	0.05-25 µg/mL	HPLC/UV	0.5 mL	None	-20	Richmond	Inactive	LC9.1
TNF-α	Human Serum	15.6-1,000 pg/mL	ELISA	0.4 mL	None	-20	Richmond	Active	ICD134
Tobramycin	Human Plasma	10-400 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P790.00
	Human Urine	100-20,000 ng/mL	LC/MS/MS	0.05 mL	None	-20	Madison	Active	P769.01
Tolbutamide	Human Plasma	0.1-100 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Richmond	Active	LCMSB226
Tolbutamide (Carboxytolbutamide) (4-Hydroxytolbutamide)	Human Urine	5-1,000 ng/mL 0.6-300 µg/mL 0.3-150 µg/mL	LC/MS/MS	0.02 mL	None	-20	Madison	Active	P738.01
	Human Plasma	0.1-100 µg/mL 5-5,000 ng/mL 2.5-2,500 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P737.02
Topiramate	Human Plasma	10-10,000 ng/mL	LC/MS/MS	0.1mL	K <sub>2</sub> EDTA	-20	Richmond	Active	LCMS470
Topiramate (See Cocktail Assay "Epilepsy")	Human Plasma	0.1-50 µg/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-70	Madison	Active	P837.00
<b>Topotecan, Total</b>	<b>Human Plasma</b>	<b>0.1-25 ng/mL</b>	<b>LC/MS/MS</b>	<b>0.05 mL</b>	<b>Sodium Heparin</b>	<b>-70</b>	<b>Richmond</b>	<b>Active</b>	<b>LCMSD507</b>

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Triamcinolone	Human Plasma	0.01-2.5 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC315
	Human Plasma	1-1,000 ng/mL	LC/MS/MS	0.1 mL	K <sub>3</sub> EDTA	-70	Richmond	Active	LCMSC315.1
Triamterene (Hydroxytriamterene Sulfate)	Human Plasma	0.5-100 ng/mL 5-1,000 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC109
	Human Urine	0.10-5 µg/mL 1-50 µg/mL	HPLC/FL	0.1 mL	None	-20	Richmond	Inactive	LC109.1
TSH	Human Serum	0.250-64 µIU/mL	RIA (IRMA)	0.4 mL	None	-20	Richmond	Inactive	ICD179
Valganciclovir (Ganciclovir)	Human Plasma	1-500 ng/mL 10-5,000 ng/mL	LC/MS/MS	0.05 mL	Sodium Citrate	-70	Richmond	Active	LCMSC410
Valproic Acid	Human Plasma	2-100 µg/mL	LC/MS/MS	0.05 mL	Sodium Heparin	-20	Madison	Active	P885.00
Venlafaxine (O-Desmethylvenlafaxine)	Human Plasma	2-200 ng/mL 2-200 ng/mL	HPLC/FL	0.5 mL	Sodium Heparin	-20	Richmond	Inactive	LC328
Verapamil (Norverapamil)	Human Plasma	1-500 ng/mL 1-500 ng/mL	HPLC/FL	1.0 mL	Sodium Heparin	-20	Madison	Inactive	P072.02
Verapamil Enantiomers (Norverapamil Enantiomers)	Human Plasma	0.5-250 ng/mL 0.5-250 ng/mL	LC/MS/MS	0.4 mL	Sodium Heparin	-20	Madison	Active	P671.00
Voriconazole (Voriconazole-N-Oxide)	Human Plasma	10-5,000 ng/mL 10-5,000 ng/mL	LC/MS/MS	0.1 mL	Sodium Heparin	-20	Madison	Active	P968.00
Warfarin Enantiomers, (unbound)	Human Plasma Ultrafiltrate	0.5-25 ng/mL 0.5-25 ng/mL	LC/MS/MS	0.05 mL	EDTA	-20	Madison	Active	P921.00

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

**Shaded entries are NEW to our list**



## LIST OF VALIDATED METHODS FEBRUARY 2010

Compound (Metabolite)	Matrix	Calibration Range	Method	Sample Volume	Anti-Coagulant	Storage Temp ° C	Lab Site	Method Status*	Method I.D.
Warfarin Enantiomers	Human Serum	5-1,000 ng/mL 5-1,000 ng/mL	LC/MS/MS	0.5 mL	None	-20	Madison	Inactive	P381.01
	Human Plasma	5-1,000 ng/mL 5-1,500 ng/mL	LC/MS/MS	0.25 mL	Sodium Heparin	-20	Richmond	Active	LCMSB194
	Human Plasma	5-1,500 ng/mL 5-1,500 ng/mL	LC/MS/MS	0.2 mL	Sodium Heparin	-20	Madison	Active	P807.01
S-warfarin (7-OH-S-warfarin)	Human Plasma	5-1,500 ng/mL 5-1,000 ng/mL	LC/MS/MS	0.2 mL	Sodium Citrate	-20	Madison	Active	P832.00
Zolpidem	Human Plasma	0.5-500 ng/mL	LC/MS/MS	0.05 mL	K <sub>2</sub> EDTA	-20	Madison	Active	P844.01
Zonisamide	Human Plasma	10-5,000 ng/mL	LC/MS/MS	0.25 mL	K <sub>3</sub> EDTA	-20	Richmond	Active	LCMS313
Zonisamide	Human Serum	10-5,000 ng/mL	LC/MS/MS	0.25 mL	None	-20	Richmond	Active	LCMS313.1

\*Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current FDA Guidance

\*\*Total = Unconjugated + Conjugated Drug.

Unless specifically noted, the free or unconjugated drug is measured.

Shaded entries are NEW to our list